

Top:

The eResponder Mobile Personal Emergency Response System (PERS) with U-TDOA Locate allows you to get help anytime, anywhere in the US with T-Mobile cellular coverage

Simple to use

1. Press the button
2. Talk to a certified emergency care specialist

Features

Works up to 2 months on a single charge
Small, lightweight and comfortable to wear
Can be worn in the shower
Works at home and away from home
Can be located with U-TDOA

Overview:

Your life should not be limited by the range of your personal emergency response system.

Your life has unlimited possibilities. With eResponder, you can feel confident staying independent and going anywhere you want to go, indoors or outdoors, because you'll know you can get help if you need it.



Go anywhere

No base station or landline phone connection is needed with eResponder. It is self-contained and has a speaker and microphone in the unit.

eResponder is your connection to help anywhere you go in the U.S. where there is T-Mobile cellular coverage ([Check coverage for your area](#)). It works both indoors and outdoors; at home and away from home; even in the shower.

You do not need to purchase any cellular carrier contracts for eResponder.



Simple to use

- 1. Press the button
- 2. Talk to an emergency care specialist
- 3. The emergency care specialist gets you the help you need.

Press one button to talk to certified emergency care specialists 24/7 if you ever need help. Emergency care specialists can send an ambulance or loved ones to you depending on what you need.

If you need to be located, emergency care specialists can connect you to 911 to utilize U-TDOA location technology.

Small, lightweight and shower-safe

eResponder weighs half as much as competitor devices.

eResponder is shower-safe so you can always keep it close at hand. However, eResponder is not designed to be submerged in water.



Long battery life

eResponder's rechargeable battery lasts up to 2 months per charge.
Other emergency alert devices available have a battery life of approximately 2 days.



eResponder vs. other personal emergency response systems

Feature	eResponder	Other Products	
		Landline	Mobile
Talk to emergency care specialists anywhere*	<input checked="" type="checkbox"/>	NO	Some do not work in shower
Emergency locate anywhere*	<input checked="" type="checkbox"/>	NO	Limitations indoors
Smallest and lightest	<input checked="" type="checkbox"/>	NO	NO
Up to 2 months per battery charge	<input checked="" type="checkbox"/>	NO	Battery lasts only 2-3 days
Easiest to use	<input checked="" type="checkbox"/>	NO	NO

* Within T-Mobile cellular coverage area in U.S.

eResponder is the smallest and lightest mobile personal emergency response system with 2-way voice communication available. It uses U-TDOA to locate users in emergencies. Unlike some mobile personal emergency response systems, eResponder can be used in the shower. However, it is not designed to be submerged in water.

How it works

How eResponder works:

1



A user needs help; the user presses the SOS button for 2 seconds

2



An emergency care specialist talks to the user via two-way voice communication

3



Help is sent to the user and the emergency contact is notified

Simple to use

1. Press the button
2. Talk to a emergency care specialist
3. The emergency care specialist gets you the help you need.

Press one button to talk to certified emergency care specialists 24/7 if you ever need help. Emergency care specialists can send an ambulance or loved ones to you depending on what you need.

If you need to be located, emergency care specialists can connect you to 911 to utilize U-TDOA location technology.

Care Specialists

What are the benefits of emergency care specialists?

Emergency care specialists are highly trained and certified and quickly respond to medical alerts 24/7.

Emergency care specialists can:

- View user-provided medical information to assist first responders in an emergency.
- Share user-provided lockbox codes with first responders in an emergency.
- Tell first responders which hospital the user prefers.


- Call caregivers, family and friends for the user.
- Help the eResponder user test their eResponder unit anytime.

Specifications

Small, lightweight and shower-safe

eResponder weighs half as much as competitor devices.

eResponder is shower-safe so you can always keep it close at hand. However, eResponder is not designed to be submerged in water.



Long battery life

eResponder's rechargeable battery lasts up to 2 months per charge.

Other mobile personal emergency response systems available have a battery life of approximately 2 days.

eResponder Specifications

Colors:	Black and white
Weight:	1.2 oz
Size:	1.75" x 1.75" x 0.69"

FAQs

1. What makes eResponder different from other personal emergency response systems?

eResponder is a very small, lightweight, and comfortable-to-wear mobile personal emergency response system (PERS) that allows users to contact emergency care specialists 24/7 from anywhere in the US with T-Mobile coverage. The small, waterproof device frees users from being tethered to their homes by speaker or landline-based systems, allowing instant communication with amazing sound quality – in the shower, indoors, outdoors- wherever the user goes.

2. How does eResponder work?

When you press the eResponder button, you are connected with an emergency care specialist. You can speak with and listen to the emergency care specialist through the built-in microphone and speaker in the eResponder device. Specialists will analyze your situation and can send EMS, loved ones or neighbors to help you, depending on what is needed.

3. Does eResponder work in the shower?

Yes, eResponder can be used in the shower. However, it is not designed to be submerged in water.

4. Why does the eResponder light turn red?

There are a number of situations in which the eResponder light turns red:

After you finish a call with an emergency care specialist: The red light flashes after five minutes to show that the eResponder is going back into standby mode. After every call, the eResponder will be powered on for five minutes just in case the emergency care specialist needs to call back the user.

When the eResponder battery is low: The red light will flash around the SOS button. When you connect the charger, the red light will flash every two seconds while charging.

If you press the SOS button and the blue light turns off, press the SOS button again: The light will turn red as eResponder dials 911.

5. Who can I contact with the eResponder?

When you press the SOS button, you will be connected with an emergency care specialist. The emergency care specialist can get you the help you need from EMS, relatives or friends.

6. How long does the eResponder battery last?

The rechargeable eResponder battery lasts up to 2 months between charging. Other emergency alert devices have a battery life of approximately 2 days.

7. How do I know when the battery is low?

When the battery is low, a red light will flash around the center button, letting you know it's time to plug it in with the included AC wall charger. The user and caregiver can also receive email and text message alerts when eResponder needs to be charged.

8. How does the eResponder help seniors and individuals with special needs?

eResponder empowers you to go anywhere in the US with T-Mobile coverage and be able to reach help. You can wear eResponder all the time, even in the shower, so you can feel confident doing your daily activities and being able to reach help if you need it.

9. How much does the eResponder unit weigh and what are the dimensions?

It weighs less than 1.2 oz. and measures 1.75"x1.75"x0.69".